



PACKAGING GUIDE


CEL/820 Basic Wastewater Laboratory


This packaging guide will help you arrange the instruments, reagents and equipment included in the CEL/820 Basic Wastewater Laboratory. The legend also serves as a materials list; use it to verify that your shipment is complete when you receive it from Hach and to re-order supplies when you need them.

Most of the chemicals and apparatus fit into the carrying case provided, except for the AccuVac® Reagent Ampuls you may want to carry with you. To transport the ampuls, a padded hinged box that holds 8 ampuls is included. Store the styrofoam shipping boxes with the remaining ampuls in a suitable place and replenish the transport box as needed.

Included Documentation:

Item Description	Cat No.
Instrument Manual, DR/800 series.....	48440-18
Procedures Manual, DR/820	48440-22





CASE BOTTOM

CEL/820 Legend (continued)

3. Demineralizer Bottle.....	14299-00	12. Ascorbic Acid Powder Pillows.....	14577-99
4. Beaker, 50-mL, low form.....	1080-41	13. DPD Tot. Chlorine Powder Pillows.....	21056-69
5. Amino Acid Reagent.....	1934-32	14. FerroVer [®] Iron Powder Pillows.....	21057-69
6. Molybdate Reagent.....	2236-32	15. NitraVer [®] 5 Powder Pillows.....	21061-69
7. Buffer Powder Pillows.....	22992-64	16. NitriVer [®] 3 Powder Pillows.....	21071-69
8. Clipper, for medium powder pillows.....	968-00	17. Thermometer, pocket, -10–110 °C.....	1877-01
9. AccuVac Transport Box		18. COD/TNT Adapter, DR/800.....	48464-00
Dissolved Oxygen AccuVac Ampuls.....	25150-25	19. DR/820 Colorimeter.....	48440-60
10. AluVer [®] Aluminum 3 Powder Pillows.....	14290-99	20. Sample Cells, 10-, 20- & 25-mL.....	24019-06
11. Bleaching 3 Powder Pillows.....	14294-49	21. Pocket Pal pH Tester.....	44350-00



FOR TECHNICAL ASSISTANCE, PRICE INFORMATION AND ORDERING:
In USA call 1-800-227-4224. Or contact the authorized HACH office or distributor serving you.